



New Blood Cell Counting Technology in the VCUHS Hematology Laboratory

Hematology Laboratory Faculty and Staff
March, 2017



On March 7, 2017, the VCUHS Clinical Laboratories will introduce a new, top-of-the line automated CBC analyzer (i.e., hematology analyzer) track line (Sysmex XN-9000) in the main Hematology Laboratory. Similar Sysmex technology will also be used in the other VCUHS laboratories providing cellular analysis, including the Emergency Department Laboratory, Stony Point Laboratory, Children's Hospital Laboratory, and Transplantation Laboratory. This change was instituted because recent technological developments in cellular analysis and computer software have led to improved analytic efficiency and new cellular analytic parameters of clinical significance. In addition to the conventional CBC parameters and five-part leukocyte differential count, the new hematology analyzers provide an extended leukocyte differential, with quantitative information regarding cells that do not normally occur in the peripheral blood, such as blasts, immature granulocytes, atypical leukocytes, and nucleated red blood cells (NRBCs). In addition, they provide additional information about reticulocytes and platelets that is of great value in patient diagnosis and management.

Sysmex Technology

The Sysmex automated CBC analyzers use the following technologies to measure and count circulating blood cells and cellular components: (1) direct current impedance, (2) advanced optical light scatter technology, (3) fluorescent flow cytometry, and (4) spectrophotometry.

The WBC, RBC, platelet counts, hemoglobin, and hematocrit are measured directly in the Sysmex instruments using a WBC channel, a RBC/platelet channel, and a separate hemoglobin channel. WBC and RBC counts are obtained using impedance technology, which is enhanced by a

sheathed stream with hydrodynamic focusing in the RBC/platelet channel, and by floating thresholds in the WBC/platelet channel to accurately discriminate cell populations. Although other CBC analyzers use the conventional modified cyanohemoglobin technique to measure hemoglobin concentration, the Sysmex instruments utilize a sodium lauryl sulfate (SLS) reagent, which forms a stable complex with oxidized hemoglobin, and is measured photometrically at 555 nm to minimize interference by turbidity and other interferences.

Various cellular indices are calculated from the direct measurements by the CBC analyzer, including the MCV, MCH, MCHC, RDW-SD, RDW-CV, MPV, and plateletcrit. The RDW-CV is the conventional RDW reported as a CV. The RDW-SD is a new parameter representing the RBC arithmetic distribution width measured at 20% of the height of the RBC curve, reported in femtoliters (fL), with a reference interval of 37 to 54 fL. The plateletcrit is the platelet volume ratio, analogous to the hematocrit. MPV is calculated from the plateletcrit and platelet count in the same manner as the erythrocyte MCV.

Fluorescent flow cytometry using a polymethine dye is used to determine the WBC, WBC differential, and the enumeration of nucleated RBCs. In the Sysmex XN-9000, the three separate channels are used for these measurements. The white cell nucleated (WNR) channel is used for WBC, NRBC, and basophil counting, while neutrophils, lymphocytes, monocytes, eosinophils, and immature granulocytes are determined from the WBC differential (WDF) channel. If blasts or abnormal lymphocytes are detected by the instrument, subsequent analysis is performed in a special white precursor cell (WPC) channel. The WBC count is automatically corrected when NRBCs are detected in the sample.




Fig. 1. Sysmex XN-9000 histograms from a patient with chronic lymphocytic leukemia (CLL). In the WDF channel (upper left) the lymphocyte population is pink, monocytes are bright green, and neutrophils are aqua. In the WNR channel (upper middle) debris is dark blue and leukocytes are aqua. RBC and PLT histograms are shown (lower left). In the RET channel, the red blood cell population is purple, and reticulocytes are pink. Mature platelets are aqua in the PLT-F channel, while the immature PLT fraction is bright green.

In addition to conventional impedance and optical methods for counting platelets, a novel fluorescent PLT (PLT-F) channel uses a fluorescent RNA dye (oxazine) in combination with extended PLT counting volume and time for platelet counting and determination of the immature platelet fraction (IPF).

Specimen Requirements, Processing, and Laboratory Analysis

The routine specimen requirement for the CBC and CBC with differential for Sysmex analysis is a 3 mL lavender top (K₂EDTA) tube. 500 uL capillary blood collection (Microtainer) tubes will be accepted for pediatric patients. The analysis will be available on all shifts and days of the week (24/7/365). Specimens should be delivered within four hours of collection. However, remote locations can be refrigerated at 20-8°C for up to 24 hours if testing cannot be completed within 4 hours of collection. Under these circumstances, transportation to the laboratory should be in a transport cooler containing cold packs.

Table 1
Adult Blood Cell Reference Ranges

Parameter	Female	Male	Units
WBC	3.9 - 11.7	3.7 - 9.7	x 10 ⁹ cells/L
RBC	3.85 - 5.16	4.54 - 5.78	x 10 ¹² cells/L
HGB	12.0 - 15.0	13.3 - 17.2	g/dL
HCT	34.8 - 45.0	38.9 - 50.9	%
MCV	78.5 - 96.4	81.2 - 94.0	fL
MCH	25.6 - 32.2	25.7 - 32.2	pg
MCHC	30.5 - 34.0	30.9 - 35.5	g/dL
RDW-CV	11.3 - 14.7	11.5 - 14.1	%
RDW-SD	36.4 - 46.3	35.1 - 43.9	fL
PLT	172 - 440	179 - 373	x 10 ⁹ cells/L
MPV	8.7 - 12.3	8.7 - 12.1	fL
IPF	0.0 - 9.9	0.0 - 9.9	%
NEU#	1.9 - 7.9	2.0-6.7	x 10 ⁹ cells/L
LYM#	1.3 - 3.6	1.1 - 3.3	x 10 ⁹ cells/L
MONO#	0.3 - 0.7	0.2 - 0.7	x 10 ⁹ cells/L
EOS#	0.0 - 0.4	0 - 0.4	x 10 ⁹ cells/L
BASO#	0.0 - 0.1	0 - 0.1	x 10 ⁹ cells/L
nRBC%	0.0 - 0.20	0.0 - 0.20	% (/100WBC)
% Retic	0.9 - 2.6	0.8 - 2.5	%
Absolute Retic Count	0.0377 - 0.1222	0.0436 - 0.1305	x 10 ⁶ cells/L
IRF	2.3 - 13.4	2.3 - 13.4	%
RET-He	32.1 - 39.1	32.1 - 39.1	pg

Continuous orders will need to be re-ordered in Cerner to continue to be performed by the new method.

Turnaround time for results is 1 hour for Stat specimens and 4 hours for routine analysis. Adult reference ranges for Sysmex blood cell counts and indices is given in the following table, and will be included with the patient result. Pediatric reference ranges will be provided with the patient report and are available on the pathology department web site at the following address (<http://www.pathology.vcu.edu/clinical-services/clinical-pathology/hematology>).

The Sysmex instruments are very sensitive to the presence of abnormal blood cells and blood cell counts generate “flags” to alert the cell counter operator to the presence of specimen abnormalities. If a significant ab-

normality is identified, a blood film is automatically prepared for further review by a technologist using digital cell morphology (Cellavision) technology and possibly a light microscope. The results may be released to the Cerner chart by the technologist, but may be referred to an attending hematopathologist for consultation if certain criteria are met. If the results are medically significant or exceed a critical value, the ordering physician will be notified by telephone and a comment entered in the Cerner result.

New FDA-Approved Reportable Cellular Parameters

The complete blood count (CBC) with differential performed in the main VCUHS Hospital Hematology Laboratory. will now include cell counts (RBC, WBC, and PLT), a five part differential count, hemoglobin, hematocrit, MCV, MCH, MCHC, RDW-CV,

reticulocyte % reticulocytes, absolute reticulocyte number, Immature Reticulocyte Fraction (IRF), Reticulocyte Hemoglobin (RET-He), Nucleated Red Blood Cell (nRBC) Count, and Immature Platelet Fraction (IPF). The IPF will not be available from specimens analyzed in the Transplantation Laboratory. The new tests are an extension of the traditional CBC and require no additional blood sample. The Immature Granulocyte (IG) count is under review by the laboratory and will be available at a later time, with notification from the laboratory.

Reticulocyte Hemoglobin (RET-He) and Immature Reticulocyte Fraction (IRF):

Reticulocyte counts are the quantity of the youngest erythrocytes normally released from the bone marrow into circulating blood. Automated analysis has led rapid and very accurate reticulocyte counting, as well as providing the reticulocyte immature reticulocyte fraction (IRF), the reticulocyte hemoglobin content (RET-He), and other parameters. The IRF assesses reticulocyte maturation by measuring the intensity of mRNA staining, with the youngest reticulocytes having the highest content. IRF has been proposed as an early marker of engraftment in bone marrow or hematopoietic stem cell transplantation and bone marrow regeneration following chemotherapy. RET-He is a reliable marker of cellular hemoglobin content that is useful in assessing the functional iron available for erythropoiesis during the previous 3–4 days. Published data show a RET-He cut-off of 29 pg/cell as an indication of iron deficiency. A value below this range is indicative of a decreased amount of iron incorporation into the RBC or iron deficiency. There is a good agreement between the RET-He and the CHr measured by other automated cell counters. These indices correlate with iron-deficient erythropoiesis and are useful markers of iron deficiency in infants and children,

adult blood donors, geriatric patients, pregnant women, and patients with chronic kidney disease undergoing hemodialysis. In addition, the RET-He is a useful in monitoring the response to iron replacement therapy and detecting iron-restricted erythropoiesis in patients receiving erythropoietin therapy. There is also intense interest in using these parameters to evaluate and monitor patients with complex anemia associated with malignancy.

Nucleated Red Blood Cells (NRBC): Nucleated red blood cells are immature erythrocytes that are commonly found in the circulation during pregnancy and the very early neonatal period. However, the finding of nRBCs at other times usually indicates markedly increased erythroid activity, or




Fig. 2. Histogram of RET channel showing location of reticulocytes in comparison to mature RBCs (blue). The three stages of reticulocyte maturation include low fluorescence reticulocytes (LFR), medium fluorescent reticulocytes (MFR), and high fluorescent reticulocytes (HFR). The combination of the two most immature stages of maturation, MFR and HFR, comprise the immature reticulocyte fraction (IRF).

Table II
Summary of Sysmex Parameters

Sysmex Parameters FDA-Approved and Reported at VCUHS		
Parameter	Technology	Units reported
Red blood cell count	Impedance with hydrodynamic focusing	RBC (x 10 ¹² cells/L)
Hemoglobin	RBC lysis by SLS, photometry at 555 nm	Hgb (g/L)
RBC indices	Impedance, calculations, and histogram analysis	Hematocrit, MCV, MCH, MCHC, RDW-CV
Reticulocyte count	Fluorescent flow cytometry	% and absolute count
Immature reticulocyte fraction	Fluorescent flow cytometry	IFR (%)
Reticulocyte hemoglobin	Fluorescent flow cytometry	RET-He (pg)
Nucleated red cell	Fluorescent flow cytometry	nRBC (%/100 WBC)
White blood cell count	Fluorescent flow cytometry	WBC (x 10 ⁹ cells/L)
WBC Five-Part Differential	Light scatter and fluorescent flow cytometry	NEU#, NEU%, LYM#, LYM%, MONO#, MONO%, EOS#, EOS%, BASO#, BASO%
Platelet count	Impedance and fluorescent flow cytometry	PLT (x 10 ⁹ cells/L)
Parameter	Technology	Units reported
Red blood cell count	Impedance with hydrodynamic focusing	RBC (x 10 ¹² cells/L)
Hemoglobin	RBC lysis by SLS, photometry at 555 nm	Hgb (g/L)
Sysmex Parameters Non-FDA Approved and/or Not Reported at VCUHS		
RDW-SD	Calculation	RDW-SD %
Immature granulocyte fraction	Fluorescent flow cytometry	IG (%)
Neutrophil granulation (NEUT-SSC)	Light scatter and fluorescent flow cytometry	NEUT-SSC
Fragmented red cells (FRC)	Fluorescent flow cytometry	FRC (# and %)
Microcytic and Macrocytic Red Blood Cell Populations	Histogram analysis	MicroR and MacroR (%)
Hypo-Hemoglobinized and Hyper-Hemoglobinized Red Blood Cells	Fluorescent flow cytometry	HYPO-He and HYPER-He (%)

damage to the bone marrow microenvironment from a hematologic neoplasm or metastatic malignancy. In these conditions, the presence of circulating nRBCs is a poor prognostic indicator. Until recently, it was virtually impossible for hematology analyzers to distinguish small mature lymphocytes from nRBCs. At that time, the only reliable method for enumerating circulating nRBCs was manual counting of peripheral blood smears, usually using a 100 cell count, with significant intra- and interobserver reproducibility and sampling error. Earlier automated nRBC enumeration by different

hematology analyzers was also limited by low sensitivity and specificity, especially when dealing with low numbers of nRBCs. In contrast, the Sysmex hematology analyzers use advanced light scattering techniques and fluorescence to rapidly and accurately detect NRBCs even at clinically significant low numbers.

Immature Platelet Fraction (IPF): Circulating immature platelets are much larger platelets that have been recently released from the bone marrow. IPF have a greater

RNA content, and can be measured by automated hematology analyzers, including the Sysmex instruments in the same manner as immature reticulocytes, and they are reported as percentage of the total platelet count (% IPF). It is well accepted that a high IPF is usually found in either consumptive or recovering thrombocytopenic disorders, while a low IPF is characteristic of bone marrow suppression disorders. Thus, the IPF is an index of thrombopoiesis and it may assist the physician in determining the cause and differential diagnosis of thrombocytopenia when used with patient information and the platelet count. IPF may help the physician determine if thrombocytopenia is due to platelet destruction or decreased platelet production. In addition, many studies have also shown that the IPF is an early indicator of marrow recovery in patients rebounding from chemotherapy or hematopoietic stem cell transplant. IPF recovery, denoted as levels >7.0%, occur on average 3.1 days earlier than platelet count recovery, and 3.8 days earlier than absolute neutrophil count recovery. Thus, the IPF may be useful to guide and possibly limit prophylactic platelet transfusions in patients undergoing marrow suppressive therapy, in view of imminent recovery of the platelet count.

Immature Granulocyte count (IG): Historically, measurement of the proportion of the circulating neutrophil population that are bands (i.e., band count) has been considered clinically useful for the diagnosis of infection, especially for neonatal sepsis. Since automated hematology analyzers could not accurately differentiate bands from neutrophilic granulocytes at other stages of maturation, manual counting using a light microscope was used for band counts. However, multi-institutional studies conducted by the College of American Pathologists and other organizations showed that band counts were inaccurate

and imprecise due to the inability of technologist to reproducibly identify band neutrophils. During the past decade, most hospitals have stopped performing manual band counts for this reason. Fortunately, other immature granulocytes (IGs) including metamyelocytes, myelocytes, and promyelocytes have better morphological definition and together can be used as an alternative to the band count. In addition, IGs are usually not detected in healthy individuals but are elevated in patients with bacterial infections, acute inflammatory disorders, cancer (marrow metastasis), tissue necrosis, acute transplant rejection, surgical and orthopedic trauma, myeloproliferative neoplasm, steroid use, and pregnancy. An increase in IGs is usually accompanied by an increase in the absolute neutrophil count (ANC), but elderly patients, neonates, and patients with myelosuppression may have elevated IGs without an elevation of the ANC. Sysmex XN hematology analyzers perform the IG count as a part of the leukocyte differential count with notably low imprecision (CV near 7%). In addition, the accuracy of these measurements compared to microscopic examination or flow cytometry with monoclonal antibodies has been shown to be high. The ability to provide a more accurate and precise automated immature granulocyte count without performing a manual differential will decrease turnaround times to provide patient results sooner. The VCUHS Hematology Laboratory is presently evaluating this parameter in our hospital.

**New Non-FDA-Approved,
Non-Reportable Cellular Parameters**

A number of measurements and calculations are performed during automated cell analysis that have not been FDA-approved, will not be reported with the patient results,

and cannot be used for *in vitro* diagnosis. These parameters should be considered experimental, and include measurements of neutrophil granularity (NEUT-SSC), fragmented red blood cells (FRC*), hypochromic red cells (%HYPO-He), hyperchromatic red cells (%HYPER-He), macrocytic red cells (MacroR), and microcytic red cells (MicroR). With appropriate IRB-approval and collaboration with the hematology laboratory faculty, these parameters can be made available for research studies in patients with anemia and other diseases. In addition, research investigations involving the RDW-SD, immature granulocyte count, immature platelet count, RET-He and other analytic parameters are greatly needed to expand knowledge in this area. The hematopathology faculty is most willing to collaborate with the faculty and housestaff of other departments with a research interest in blood cell counts and parameters.

Neutrophil Granulation (NEUT-SSC):

Light scattered at a 90° angle from the axis of a laser beam (side scatter, SSC) is altered by the internal complexity of cells passing through the laser beam and analysis of SSC provides information about the cytoplasmic density and number of cytoplasmic granules in a cell population. Among the peripheral blood leukocyte population, neutrophils and eosinophils exhibit the highest SSC, in comparison the monocytes and lymphocytes. The Sysmex cell counters analyze the side-scattered light (SSC) of the WDF channel to generate the parameter NEUT-SSC, which is a measure of the granularity of the neutrophil population. A low NEUT-SSC value indicates a hypogranular neutrophil population, while the NEUT-SSC value is increased in hypergranularity. Neutrophil hypogranularity is a common feature of neutrophil dysplasia found in the myelodysplastic syndromes (MDS) and some myeloproliferative disor-

ders such as chronic myelomonocytic leukemia (CMML) and atypical chronic myeloid leukaemia (aCML). Neutrophil hypergranularity is exhibited in the reactive neutrophils found in patients with an infection or acute phase reaction. The automated detection of hypogranular neutrophils is potentially of great diagnostic value in differentiating between MDS and reactive and benign idiopathic and hereditary causes of neutrophilia, especially when used in conjunction with absolute cell counts, the immature granulocyte count, and other parameters.

Fragmented Red Blood Cells (FRC*):

The Sysmex cell counters analyze a specific area of the reticulocyte (RET) histogram to determine the number and percentage of fragmented red blood cells (FRC% and FRC#). Fragmented red blood cells are usually a consequence of mechanical damage induced by turbulent blood flow or contact of the red cells with a pathologically altered endothelium. These abnormal shear forces damage the red blood cells, producing cell remnants that appear as ‘helmet’ cells, schistocytes and other odd shapes that are collectively termed “fragmented red cells” when viewed under a microscope. The finding of fragmented red cells by manual light microscopy has been used for decades to differentiate disseminated intravascular coagulation (DIC), thrombotic thrombocytopenia purpura, and similar diseases from other causes of hemolytic anemia. The subjective automated measurement of this parameter should greatly increase the diagnostic efficacy of hemolytic anemia diagnosis. Forthcoming advanced red cell morphology analytic software for the Cellvision digital morphology analyzer should complement the Sysmex FRC parameter.

Microcytic and Macrocytic Red Blood

Cell Populations: Measurement of the proportion of the red blood cell population that is microcytic and macrocytic is determined in the Sysmex cell analyzers by analysis of the RBC size distribution (i.e., MCV) histogram produced by the RBC/PLT channel. In this channel, hydrodynamic focusing with sheath flow direct current) impedance is used to count RBCs and platelets (PLT) and produce the histogram. Normally, this histogram shows a nearly Gaussian distribution, but may be skewed in patients with microcytosis or macrocytosis. The percentage of microcytic (MicroR) and macrocytic (MacroR) cells is determined by placing discriminators at the lower and upper area of the histogram. These parameters are helpful in refining the diagnosis of anemia, since the MCV can be normal, even in the presence of dimorphic red cell populations or increased microcytic or macrocytic red cells.

Hypo-Hemoglobinized (HYPO-He) and Hyper-Hemoglobinized (HYPER-He) Red

Blood Cells: HYPO-He is the percentage of RBC with a cellular hemoglobin content lower than 17 pg, whereas HYPER-He is the percentage of RBC with a cellular hemoglobin content higher than 49 pg. These parameters are determined in the reticulocyte (RET) channel, where they are derived from a measurement of the hemoglobin content of all mature RBC (RBC-H_e). In conjunction with the reticulocyte count, IRF, RET-He and other parameters, the HYPO-He and HYPER-He are useful in predicting, diagnosing, and monitoring iron deficiency in patients with chronic renal failure and other diseases.

Other Parameters: The Sysmex XN-9000 analyzer is capable of providing an extended leukocyte differential count with as many as 22 parameters that are generated with the CBC. In addition to the parameters discussed above, these include neutrophil

complexity and width of dispersion (NE-WX), neutrophil fluorescence intensity (NE-SFL), neutrophil fluorescence intensity and the width of dispersion (NE-WY), neutrophil cell size (NE-FSC), neutrophil cell size and the width of dispersion (NE-WZ), lymphocyte cell complexity (LY-X), lymphocyte complexity and width of dispersion (LY-WX), lymphocyte fluorescence intensity (LY-Y), lymphocyte fluorescence intensity and the width of dispersion (LY-WY), lymphocyte cell size (LY-Z), lymphocyte cell size and the width of dispersion (LY-WZ), monocyte cell complexity (MO-X), monocyte complexity and width of dispersion (MO-WX), monocyte fluorescence intensity (MO-Y), monocyte fluorescence intensity and the width of dispersion (MO-WY), monocyte cell size (MO-Z), and monocyte cell size and the width of dispersion (MO-WZ).

Additional information

Copies of this brochure are available on the VCUHS hematology laboratory web site (<http://www.pathology.vcu.edu/clinical-services/clinical-pathology/hematology>). Additional information can be obtained from Client Services at 804-828-PATH, or from the hematopathology faculty. A complete bibliography in EndNote format is available from Dr. Riley (roger.riley@vcuhealth.org). The following webinars are also available:

Measurement of immature cells in peripheral blood (<https://www.youtube.com/watch?v=z7IGJqSLuAw>)
Utility of RET-He (<https://www.youtube.com/watch?v=93gYGZ7CWEY>)
Utility of IPF (<https://www.youtube.com/watch?v=zcBxQnDDtAc>)
Utility of IG (<https://www.youtube.com/watch?v=bBG1O-0zca0>)
Case studies demonstrating the clinical application of the advanced clinical parameters (<https://www.sysmex.com/us/en/Education/Webinars/Pages/LearnMore.aspx?WebinarID=67&Upcoming=0>)

Acute anemia toolkit: bloodless strategies to optimize patient outcomes
<https://www.sysmex.com/us/en/Education/Webinars/Pages/LearnMore.aspx?WebinarID=74&Upcoming=0>
A clinical overview and discussion of laboratory medicine in an oncology practice
<http://www.sysmexamerica.com/CRC/Webinars/Pages/LearnMore.aspx?WebinarID=63&Upcoming=0>

Selected References

Interpretation of Blood Cell Counts

- Bain B.J. Diagnosis from the blood smear. *N. Engl. J. Med.* 353(5):498-507, 2005.
- Briggs C. et al. ICSH Guideline for worldwide point-of-care testing in haematology with special reference to the complete blood count. *Int. J. Lab. Hematol.* 30(2):105-116, 2008.
- Gulati G., Song J., Florea A.D., Gong J. Purpose and criteria for blood smear scan, blood smear examination, and blood smear review. *Ann. Lab. Med.* 33(1):1-7, 2013.
- Jakubik L.D. et al. The ABCs of pediatric laboratory interpretation: understanding the CBC with differential and LFTs. *Pediatr. Nurs.* 29(2):97-103, 2003.
- Leach M. Interpretation of the full blood count in systemic disease--a guide for the physician. *J. R. Coll. Physicians Edinb.* 44(1):36-41, 2014.
- Milcic T.L. The complete blood count. *Neonatal Netw.* 29(2):109-115, 2010.
- Tefferi A., Hanson C.A., Inwards D.J. How to interpret and pursue an abnormal complete blood cell count in adults. *Mayo Clin. Proc.* 80(7):923-936, 2005.
- Thachil J. Do we need 'routine' blood counts? *Br. J. Hosp. Med. (Lond).* 75(11):644-646, 2014.
- Valent P. Low blood counts: immune mediated, idiopathic, or myelodysplasia. *Hematology Am. Soc. Hematol. Educ. Program.* 2012(485-491), 2012.

Cell Counting Technology

- Becker P.H. et al. Performance evaluation of the Sysmex XN-1000 hematology analyzer in assessment of the white blood cell count differential in pediatric specimens. *Int. J. Lab. Hematol.* 38(1):54-63, 2016.
- Buttarelo M. and Plebani M. Automated blood cell counts: state of the art. *Am. J. Clin. Pathol.* 130(1):104-116, 2008.
- Eilertsen H. and Hagve T.A. Do the flags related to immature granulocytes reported by the Sysmex XE-5000 warrant a microscopic slide review? *Am. J. Clin. Pathol.* 142(4):553-560, 2014.
- Furundarena J.R. et al. Comparison of abnormal cell flagging of the hematology analyzers Sysmex XN and Sysmex XE-5000 in oncohematologic patients. *Int. J. Lab. Hematol.* 39(1):58-67, 2017.
- Hagner R. The manual differential enters the digital age. *MLO Med Lab Obs.* 44(5):20-21, 2012.
- Hotton J. et al. Performance and abnormal cell flagging comparisons of three automated blood cell counters: Cell-Dyn Sapphire, DxH-800, and XN-2000. *Am. J. Clin. Pathol.* 140(6):845-852, 2013.
- International Council for Standardization in Haematology W.G., Briggs C., Culp N., et al. ICSH guidelines for the evaluation of blood cell analyzers including those used for differential leucocyte and reticulocyte counting. *Int. J. Lab. Hematol.* 36(6):613-627, 2014.
- Lecompte, T.P. and Bernimoulii, M.P. Novel parameters in blood cell counters. *Clin Lab Med.* 35(1):209-224, 2015.
- Meintker L. et al. Comparison of automated differential blood cell counts from Abbott Sapphire, Siemens Advia 120, Beckman Coulter DxH 800, and Sysmex XE-2100 in normal and pathologic samples. *Am. J. Clin. Pathol.* 139(5):641-650, 2013.
- Noordegraaf M. et al. Diagnostic efficiency of the Sysmex XN WPC channel for the reduction of blood smears. *Clin. Biochem.* 49(16-17):1292-1294, 2016.

Seo J.Y. et al. Performance evaluation of the new hematology analyzer Sysmex XN-series. *Int J Lab Hematol.* 37(2):155-164, 2015.

Anemia and Reticulocyte Analysis

Ageeli, A.A. et al. Reticulocyte hemoglobin content and iron deficiency: a retrospective study in adults. *Genet. Test Mol. Biomarkers.* 17(4):278-283, 2013.

Balci, Y.I. et al. Evaluation of Reticulocyte Parameters in Iron Deficiency, Vitamin B12 Deficiency and Mixed Anemia. *Clin. Lab.* 62(3):343-347, 2016.

Buttarello, M. et al. Evaluation of the hypochromic erythrocyte and reticulocyte hemoglobin content provided by the Sysmex XE-5000 analyzer in diagnosis of iron deficiency erythropoiesis. *Clin. Chem. Lab. Med.* 54(12):1939-1945, 2016.

Darveau, M. and Lachance, P. Reticulocyte hemoglobin content in critically ill patients. *Anesthesiology.* 113(6):1479, 2010.

Davidkova, S. et al. Comparison of reticulocyte hemoglobin equivalent with traditional markers of iron and erythropoiesis in pediatric dialysis. *Pediatr. Nephrol.* 31(5):819-826, 2016.

International Council for Standardization in Haematology W.G., Briggs C., Culp N., et al. ICSH guidelines for the evaluation of blood cell analyzers including those used for differential leucocyte and reticulocyte counting. *Int. J. Lab. Hematol.* 36(6):613-627, 2014.

Lesesve, J.F. et al. Fragmented red blood cells automated measurement is a useful parameter to exclude schistocytes on the blood film. *Int. J. Lab. Hematol.* 34(6):566-576, 2012.

Lesesve, J.F. et al. Fragmented red cells reference range for the Sysmex XN(R)-series of automated blood cell counters. *Int. J. Lab. Hematol.* 37(5):583-587, 2015.

Lopez-Ruzafa, E. et al. Reference Values of Reticulocyte Hemoglobin Content and Their Relation With Other Indicators of Iron Status in Healthy Children. *J. Pediatr. Hematol. Oncol.* 38(7):e207-212, 2016.

Lorenz, L. et al. Reticulocyte haemoglobin content as a marker of iron deficiency. *Arch. Dis.*

Child Fetal Neonatal Ed. 100(3):F198-202, 2015.

Mast, A.E. et al. Reticulocyte hemoglobin content. *Am. J. Hematol.* 83(4):307-310, 2008.

Mateos, M.E. et al. Reticulocyte hemoglobin content for the diagnosis of iron deficiency. *J. Pediatr. Hematol. Oncol.* 30(7):539-542, 2008.

Meier, E.R. et al. Absolute Reticulocyte count acts as a surrogate for fetal hemoglobin in infants and children with sickle cell anemia. *PLoS One.* 10(9):e0136672, 2015.

Miwa, N. et al. Usefulness of measuring reticulocyte hemoglobin equivalent in the management of haemodialysis patients with iron deficiency. *Int. J. Lab. Hematol.* 32(2):248-255, 2010.

Morkis, I.V. et al. Assessment of immature platelet fraction and immature reticulocyte fraction as predictors of engraftment after hematopoietic stem cell transplantation. *Int. J. Lab. Hematol.* 37(2):259-264, 2015.

Osta, V. et al. Utility of new mature erythrocyte and reticulocyte indices in screening for iron-deficiency anemia in a pediatric population. *Int. J. Lab. Hematol.* 35(4):400-405, 2013.

Park, S.H. et al. Evaluation of parameters obtained from the Sysmex XN-2000 for predicting the recovery of the absolute neutrophil count and platelets after hematopoietic stem cell transplantation. *Int. J. Lab. Hematol.* 38(2):198-208, 2016.

Parodi, E. et al. Absolute reticulocyte count and reticulocyte hemoglobin content as predictors of early response to exclusive oral iron in children with iron deficiency anemia. *Anemia.* 2016(7345835), 2016.

Piva, E. et al. Automated reticulocyte counting: state of the art and clinical applications in the evaluation of erythropoiesis. *Clin. Chem. Lab. Med.* 48(10):1369-1380, 2010.

Piva, E. et al. Clinical utility of reticulocyte parameters. *Clin. Lab. Med.* 35(1):133-163, 2015.

Piva, E. et al. Automated reticulocyte counting: state of the art and clinical applications in the evaluation of erythropoiesis. *Clin. Chem. Lab. Med.* 48(10):1369-1380, 2010.

Piva, E. et al. Clinical utility of reticulocyte parameters. *Clin. Lab. Med.* 35(1):133-163, 2015.

Priya, P.P. et al. Role of absolute reticulocyte count in evaluation of pancytopenia-a hospital based study. *J. Clin. Diagn. Res.* 8(8):FC01-03, 2014.

Proto-Siqueira, R. Reticulocyte parameters: why should clinical laboratories evaluate and report them? *Rev Bras Hematol Hemoter.* 36(1):3-4, 2014.

Raja-Sabudin, R.Z. et al. Immature reticulocyte fraction is an early predictor of bone marrow recovery post chemotherapy in patients with acute leukemia. *Saudi Med. J.* 35(4):346-349, 2014.

Rauf, S.E. et al. Immature reticulocyte fraction and absolute neutrophil count as predictor of hemopoietic recovery in patients with acute lymphoblastic leukemia on remission induction chemotherapy. *Turk. J. Haematol.* 33(2):131-134, 2016.

Riley, R.S. et al. Reticulocyte analysis by flow cytometry and other techniques. *Hematol. Oncol. Clin. North Am.* 16(2):373-420, 2002.

Riley, R.S. et al. Reticulocytes and reticulocyte enumeration. *J. Clin. Lab. Anal.* 15(5):267-294, 2001.

Schapkaitz, E. et al. Diagnosis of iron deficiency anaemia in hospital patients: Use of the reticulocyte haemoglobin content to differentiate iron deficiency anaemia from anaemia of chronic disease. *S. Afr. Med. J.* 106(1):53-54, 2015.

Schapkaitz, E. et al. Diagnosis of iron deficiency anaemia in hospital patients: Use of the reticulocyte haemoglobin content to differentiate iron deficiency anaemia from anaemia of chronic disease. *S. Afr. Med. J.* 106(1):53-54, 2015.

Schoorl, M. et al. Efficacy of advanced discriminating algorithms for screening on iron-deficiency anemia and beta-thalassemia trait: a multicenter evaluation. *Am. J. Clin. Pathol.* 138(2):300-304, 2012.

Thomas, L. et al. Reticulocyte hemoglobin measurement--comparison of two methods in the diagnosis of iron-restricted erythropoiesis. *Clin. Chem. Lab. Med.* 43(11):1193-1202, 2005.

Urrechaga, E. et al. Erythrocyte and reticulocyte indices in the assessment of erythropoiesis activity and iron availability. *Int. J. Lab. Hematol.* 35(2):144-149, 2013.

Urrechaga, E. et al. Percentage of hypochromic erythrocytes and reticulocyte hemoglobin equivalent predictors of response to intravenous iron in hemodialysis patients. *Int. J. Lab. Hematol.* 38(4):360-365, 2016.

Urrechaga E. et al. Erythrocyte and reticulocyte parameters in iron deficiency and thalassemia. *J. Clin. Lab. Anal.* 25(3):223-228, 2011.

Urrechaga, E. et al. Percentage of hypochromic erythrocytes and reticulocyte hemoglobin equivalent predictors of response to intravenous iron in hemodialysis patients. *Int J Lab Hematol.* 38(4):360-365, 2016.

Weimann, A. et al. Delta-He, Ret-He and a new diagnostic plot for differential diagnosis and therapy monitoring of patients suffering from various disease-specific types of anemia. *Clin. Lab.* 62(4):667-677, 2016.

Wollmann, M. et al. Reticulocyte maturity indices in iron deficiency anemia. *Rev Bras Hematol Hemoter.* 36(1):25-28, 2014.

White Blood Cells and Immature Granulocytes

Ansari-Lari, M.A. et al. Immature granulocyte measurement using the Sysmex XE-2100. Relationship to infection and sepsis. *Am. J. Clin. Pathol.* 120(5):795-799, 2003.

Arneth, B.M et al. Novel parameters of extended complete blood cell count under fluorescence flow cytometry in patients with sepsis. *J. Clin. Lab. Anal.* 28(2):130-135, 2014.

Bartels, P.C. et al. Evaluation of Sysmex SF-3000 performance concerning interpretive morphology flagging of the leucocyte differential count. *Clin. Lab. Haematol.* 19(3):187-190, 1997.

Bernstein, L.H. and Rucinski, J. Measurement of granulocyte maturation may improve the early diagnosis of the septic state. *Clin. Chem. Lab. Med.* 49(12):2089-2095, 2011.

Briggs, C. et al. New quantitative parameters on a recently introduced automated blood cell counter--the XE 2100. *Clin. Lab. Haematol.* 22(6):345-350, 2000.

Buoro, S. et al. Extended leukocyte differential count and C-reactive protein in septic patients with liver impairment: diagnostic approach to evaluate sepsis in intensive care unit. *Ann. Transl. Med.* 3(17):244, 2015.

Buttarello, M. et al. Sysmex SE-9000 hematology analyzer: performance evaluation on leukocyte differential counts using an NCCLS H20-A protocol. *National Committee for Clinical Laboratory Standards. Am. J. Clin. Pathol.* 108(6):674-686, 1997.

Chabot-Richards, D.S. and George, T.I. White blood cell counts: reference methodology. *Clin. Lab. Med.* 35(1):11-24, 2015.

Che, Y.Q. et al. Identification of immature granulocytes in cancer chemotherapy patients by cell counting vs. microscopic examination of blood smears. *Mol. Clin. Oncol.* 2(2):207-211, 2014.

Cimenti, C. et al. The predictive value of immature granulocyte count and immature myeloid information in the diagnosis of neonatal sepsis. *Clin. Chem. Lab. Med.* 50(8):1429-1432, 2012.

Cornet, E. et al. Contribution of the new XN-1000 parameters NEUT-RI and NEUT-WY for managing patients with immature granulocytes. *Int. J. Lab. Hematol.* 37(5):e123-126, 2015.

Eilertsen, H. and Hagve, T.A. Do the flags related to immature granulocytes reported by the Sysmex XE-5000 warrant a microscopic slide review? *Am. J. Clin. Pathol.* 142(4):553-560, 2014.

Fernandes, B. and Hamaguchi, Y. Automated enumeration of immature granulocytes. *Am. J. Clin. Pathol.* 128(3):454-463, 2007.

Field, D. et al. Performance evaluation of the immature granulocyte parameter on the Sysmex XE-2100 automated hematology analyzer. *Lab. Hematol.* 12(1):11-14, 2006.

Ha, S.O. et al. Fraction of immature granulocytes reflects severity but not mortality in sepsis. *Scand. J. Clin. Lab. Invest.* 75(1):36-43, 2015.

Maenhout, T.M. and Marcelis, L. Immature granulocyte count in peripheral blood by the Sysmex

haematology XN series compared to microscopic differentiation. *J. Clin. Pathol.* 67(7):648-650, 2014.

Meintker, L. et al. Comparison of automated differential blood cell counts from Abbott Sapphire, Siemens Advia 120, Beckman Coulter DxH 800, and Sysmex XE-2100 in normal and pathologic samples. *Am J Clin Pathol.* 139(5):641-650, 2013.

Nierhaus, A. et al. Revisiting the white blood cell count: immature granulocytes count as a diagnostic marker to discriminate between SIRS and sepsis--a prospective, observational study. *BMC Immunol.* 14(8), 2013.

Nigro, K.G. et al. Performance of an automated immature granulocyte count as a predictor of neonatal sepsis. *Am. J. Clin. Pathol.* 123(4):618-624, 2005.

Park, S.H. et al. Evaluation of parameters obtained from the Sysmex XN-2000 for predicting the recovery of the absolute neutrophil count and platelets after hematopoietic stem cell transplantation. *Int. J. Lab. Hematol.* 38(2):198-208, 2016.

Roehrl, M.H. and Wang, J.Y. Immature granulocytes in pregnancy: a story of Virchow, anxious fathers, and expectant mothers. *Am. J. Hematol.* 86(3):307-308, 2011.

Senthilnayagam, B. et al. Automated measurement of immature granulocytes: performance characteristics and utility in routine clinical practice. *Patholog. Res. Int.* 2012(483670), 2012.

van der Geest, P.J. et al. Immature granulocytes predict microbial infection and its adverse sequelae in the intensive care unit. *J. Crit. Care.* 29(4):523-527, 2014.

Platelets and IPF

Abe, Y. et al. A simple technique to determine thrombopoiesis level using immature platelet fraction (IPF). *Thromb. Res.* 118:463-469, 2006.

Adly, A.A. et al. Evaluation of the immature platelet fraction in the diagnosis and prognosis of childhood immune thrombocytopenia. *Platelets* 26(7):645, 2015.

Arneth, B. et al. Technology and new fluorescence flow cytometry parameters in hematological analyzers. *J. Clin. Lab. Anal.* 29(3):175, 2015.

Barsam, S.J. et al. Platelet production and platelet destruction: assessing mechanisms of treatment effect in immune thrombocytopenia. *Blood* 117:5723–5732, 2011.

Bat, T. et al. Measurement of the absolute immature platelet number reflects marrow production and is not impacted by platelet transfusion. *Transfusion* 53(6):1201, 2013.

Briggs, C. et al. Assessment of an immature platelet fraction (IPF) in peripheral thrombocytopenia. *Br. J. Haematol.* 126:93–99, 2004.

Briggs, C. et al. Immature platelet fraction measurement: a future guide to platelet transfusion requirement after haematopoietic stem cell transplantation. *Transfus. Med.* 16:101–109, 2006.

Briggs C. et al. Immature platelet fraction measurement: a future guide to platelet transfusion requirement after haematopoietic stem cell transplantation. *Transfus. Med.* 16(2):101-109, 2006.

Buoro, S. et al. Abnormal leukocyte scattergrams and immature platelet fraction on Sysmex XN-9000 analyzer: a new diagnostic tool for altered megakaryopoiesis? *Scand. J. Clin. Lab. Invest.* 77(1):73-75, 2017.

Cesari, F. et al. High platelet turnover and reactivity in renal transplant recipients patients. *Thrombosis and Haemostasis* 104:804–810, 2010.

Cremer, M. et al. Immature platelet fraction as novel laboratory parameter predicting the course of neonatal thrombocytopenia. *Br. J. Haematol.* 144:619–621, 2009.

Di Mario, A. et al. Immature platelet fraction (IPF) in hospitalized patients with neutrophilia and suspected bacterial infection. *J. Infect.* 59(3):201-206, 2009.

Everett, T.R. et al. Immature platelet fraction analysis demonstrates a difference in thrombopoiesis between normotensive and preeclamptic

pregnancies. *Thromb Haemost* 111(6):1177, 2014.

Geara, C. et al. Comparative study of quantitative performances between the new Sysmex XN-L (XN-550) haematology analyser and the XN-9000 in a routine laboratory. *Int. J. Lab. Haematol.* 38(1):e10, 2016.

Goncalo, A. et al. Predictive value of immature reticulocyte and platelet fractions in hematopoietic recovery of allograft patients. *Transplant Proc* 43:241–243, 2011.

Greene, L.A. et al. Beyond the platelet count: immature platelet fraction and thromboelastometry correlate with bleeding in patients with immune thrombocytopenia. *Br J Haematol*; 166(4):592, 2015.

Hoffmann, J.J. Reticulated platelets: analytical aspects and clinical utility. *Clin. Chem. Lab. Med.* 52(8):1107-1117, 2014.

Ibrahim, H. et al. (2014): Association of immature platelets with adverse cardiovascular outcomes. *J Am Coll Cardiol* 64:2122, 2014.

Jaing, T.H. et al. Assessment of platelet activation and immature platelet fraction as predictors of platelet engraftment after hematopoietic stem cell transplantation. *Cell Transplant.* 25:1259, 2016.

Kickler, T.S. et al. A clinical evaluation of high fluorescent platelet fraction percentage in thrombocytopenia. *Am. J. Clin. Pathol.* 125(2):282-287, 2006.

Ko, Y. et al. Establishment of reference interval for immature platelet fraction. *Int. J. Lab. Haematol.* 35(5): 528-533, 2013.

Lecompte, T.P. and Bernimoulin, M.P. Novel parameters in blood cell counters. *Clin. Lab. Med.* 35(1):209-224, 2015.

Mao, W. et al. Immature platelet fraction values predict recovery of platelet counts following liver transplantation. *Clin. Res. Hepatol. Gastroenterol.* 39(4):469, 2015.

Miyazaki, K. et al. Immature platelet fraction measurement is influenced by platelet size and is a useful parameter for discrimination of macrothrombocytopenia. *Hematology.* 20(10):587-592, 2015.

Moraes, D. et al. Immature platelet fraction in hypertensive pregnancy. *Platelets*. 27(4):333, 2016.

Morkis, I.V.C. et al. Assessment of immature platelet fraction and immature reticulocyte fraction as predictors of engraftment after hematopoietic stem cell transplantation. *Int J Lab Hematol* 37(2): 259, 2015.

Naz, A. et al. Importance of immature platelet fraction as predictor of immune thrombocytopenic purpura. *Pak. J. Med. Sci.* 32(3):575-579, 2016.

Park, S.H. et al. Evaluation of parameters obtained from the Sysmex XN-2000 for predicting the recovery of the absolute neutrophil count and platelets after hematopoietic stem cell transplantation. *Int. J. Lab. Hematol.* 38(2):198-208, 2016.

Park, SH et al. (2014): The Sysmex XN-2000 hematology autoanalyzer provides a highly accurate platelet count than the former Sysmex XE-2100 system based on comparison with the CD41/CD61 immunoplatelet reference method of flow cytometry. *Ann Lab Med.* 34(6):471-474, 2014.

Saigo K. et al. Usefulness of immature platelet fraction for the clinical evaluation of myelodysplastic syndromes. *Lab. Hematol.* 15(2):13-16, 2009.

Sakuragi, M. et al. Clinical significance of IPF% or RP% measurement in distinguishing primary immune thrombocytopenia from aplastic thrombocytopenic disorders. *Int. J. Hematol.* 101(4):369, 2015.

Schoorl, M et al. New fluorescent method (PLT-F) on Sysmex XN2000 hematology analyzer achieved higher accuracy in low platelet counting. *Am. J. Clin. Pathol.* 140:495-499, 2013.

Seo, J.Y. et al. Performance evaluation of the new hematology analyzer Sysmex XN-series. *Int J Lab Hematol.* 37(2):155, 2015.

Sinclair, L. The immature platelet fraction: an assessment of its application to a routine clinical laboratory. *Aust. J. Med. Sci.* 33(2):48, 2012.

Strauss, G. et al. Immature Platelet Count: a simple parameter for distinguishing thrombocytopenia in pediatric acute lymphocytic leukemia from immune thrombocytopenia. *Pediatr Blood Cancer* 57(4):641-647, 2010.

Takami, A. et al. Immature platelet fraction for prediction of platelet engraftment after allogeneic stem cell transplantation. *Bone Marrow Transplant* 39:501-507, 2007.

Tanaka, Y. et al. Performance evaluation of platelet counting by novel fluorescent dye staining in the XN-series automated hematology analyzers. *J. Clin. Lab. Anal.* 28(5):341, 2014.

Van der Linden, N. et al. Immature platelet fraction (IPF) measured on the Sysmex XN haemocytometer predicts thrombopoietic recovery after autologous stem cell transplantation. *Eur J. Haematol.* 93(2):150, 2014.

Zucker ML et al. Mechanism of thrombocytopenia in chronic hepatitis C as evaluated by the immature platelet fraction. *Int. J. Lab. Hematol.* 34:525-532, 2012.

Body fluid analysis

Boer, K. et al. Evaluation of the XE-5000 for the automated analysis of blood cells in cerebrospinal fluid. *Clin. Biochem.* 42(7-8):684-691, 2009.

Cho, Y.U. et al. Body fluid cellular analysis using the Sysmex XN-2000 automatic hematology analyzer: focusing on malignant samples. *Int. J. Lab. Hematol.* 37(3):346-356, 2015.

De Jonge, R. et al. Automated analysis of pleural fluid total and differential leukocyte counts with the Sysmex XE-2100. *Clin. Chem. Med. Lab.* 44: 1367-1371, 2006.

De Jonge, R. et al. Evaluation of the new body fluid mode on the Sysmex XE-5000 for counting leukocytes and erythrocytes in cerebrospinal fluid and other body fluids. *Clin. Chem. Lab. Med.* 48: 665-675, 2010.

De Jonge R et al. Evaluation of the new body fluid mode on the Sysmex XE-5000 for counting leukocytes and erythrocytes in cerebrospinal fluid and other body fluids. *Clin. Chem. Lab. Med.* 48: 665-675, 2010.

Fleming, C. et al. Validation of the body fluid module on the new Sysmex XN-1000 for counting blood cells in cerebrospinal fluid and other body fluids. *Clin. Chem. Med. Lab.* 50: 1791-1798, 2012.

Fleming, C. et al. Improved software on the Sysmex XE-5000 BF mode for counting leukocytes in cerebrospinal fluid. *Clin Chem Lab Med* 51: e61-63, 2013.

Fleming, C. et al. Clinical relevance and contemporary methods for counting blood cells in body fluids suspected of inflammatory disease. *Clin. Chem. Lab. Med.* 53(11):1689,, 2015.

Genc, S. et al. Evaluation of cell counting in body fluids: comparison of two automated hematology analyzers with manual microscopy. *Clin. Lab.* 62(12):2449-2453, 2016.

Labaere, D. et al. Detection of malignant cells in serous body fluids by counting high-fluorescent cells on the Sysmex XN-2000 hematology analyzer. *Int. J. Lab. Hematol.* 37(5):715, 2015.

Labaere, D. et al. Detection of malignant cells in serous body fluids by counting high-fluorescent cells on the Sysmex XN-2000 hematology analyzer. *Int. J. Lab. Hematol.* 37(5):715-722, 2015.

Li, A. et al. Automated white blood cell counts in cerebrospinal fluid using the body fluid mode on the platform Sysmex XE-5000. *Scand. J. Clin. Lab. Invest.* 74(8): 673, 2014.

Paris, A. et al. Performance evaluation of the body fluid mode on the platform Sysmex XE-5000 series automated hematology analyzer. *Int. J. Lab. Hematol.* 32: 539- 547, 2010.

Paris, A. et al. Performance evaluation of the body fluid mode on the platform Sysmex XE-5000 series automated hematology analyzer. *Int J Lab Hematol* 32: 539- 547, 2010.

Riedl, J.A. et al. Automated morphological analysis of cells in body fluids by the digital microscopy system DM96. *J. Clin. Pathol.* 63: 538-543, 2010.

Sandhaus, L.M. Body fluid cell counts by automated methods. *Clin. Lab. Med.* 35(1):93-103, 2015.

Sandhaus, L.M. et al. Automated cerebrospinal fluid cell counts using the Sysmex XE-5000: is it

time for new reference ranges? *Am. J. Clin. Pathol.* 134: 734-738, 2010.

Seo, J.Y. et al. Performance evaluation of the new hematology analyzer Sysmex XN-series. *Int. J. Lab. Hematol.* 37(2):155-164, 2015.

Tanaka, Y. et al. Evaluation of the body fluid mode of automated hematology analyzer XN-series for extremely low peripheral white blood cell counts. *Int. J. Lab. Hematol.* 36(1):e3-7, 2014.

Zimmermann, M. et al. Automated vs. manual cerebrospinal fluid cell counts: a work and cost analysis comparing the Sysmex XE-5000 and the Fuchs-Rosenthal manual counting chamber. *Int. J. Lab. Hematol.* 33: 629-637, 2011.

NRBCs

Bruegel M., Nagel D., Funk M., Fuhrmann P., Zander J., Teupser D. Comparison of five automated hematology analyzers in a university hospital setting: Abbott Cell-Dyn Sapphire, Beckman Coulter DxH 800, Siemens Advia 2120i, Sysmex XE-5000, and Sysmex XN-2000. *Clin Chem. Lab. Med.* 53(7):1057-1071, 2015.

Danise P., Maconi M., Barrella F., et al. Evaluation of nucleated red blood cells in the peripheral blood of hematological diseases. *Clin. Chem. Lab. Med.* 50(2):357-360, 2011.

Hwang D.H., Dorfman D.M., Hwang D.G., Senna P., Pozdnyakova O. Automated nucleated RBC measurement using the Sysmex XE-5000 hematology analyzer: Frequency and Clinical Significance of the Nucleated RBCs. *Am. J. Clin. Pathol.* 145(3):379-384, 2016.

Tantanate C., Klinbua C. Performance evaluation of the automated nucleated red blood cell enumeration on Sysmex XN analyser. *Int. J. Lab. Hematol.* 37(3):341-345, 2015.

Wang F.S., Itose Y., Tsuji T., Hamaguchi Y., Hirai K., Sakata T. Development and clinical application of nucleated red blood cell counting and staging on the automated haematology analyser XE-2100. *Clin. Lab. Haematol.* 25(1):17-23, 2003.